Hersteller

lege artis Pharma GmbH + Co. KG
Postfach 60, D-72132 Dettenhausen
Breitwasenring 1, D-72135 Dettenhausen
Telefon +49 (0) 71 57 / 56 45 - 0
Fax +49 (0) 71 57/ 56 45 50
E-Mail: info@legeartis.de
www.legeartis.de

Bezeichnung

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Name

HISTOLITH NaOCl 1%
For rinsing of root canals
For dental use only

Composition

1 ml of solution contains purified water, 10.5 mg (1.05% m/V) sodium hypochlorite (corresponds to 1.0% m/V active chlorine), sodium chloride and sodium hydroxide

Indications

For rinsing and cleaning of root canals during root canal preparation.

Contra-indications

Allergies to chlorine
Open apical foramen

Side-effects

When used as intended in the root canal, no side effects are known.
If sodium hypochlorite passes through the apex, irritations of the periapical tissue are possible. On healthy tissue, HISTOLITH NaOCl 1% has an irritant effect.
In the event of serious incidents in connection with the product, please contact the manufacturer and the responsible authority immediately.

Notices and Precautions

Caution, irritant.
Avoid contact of sodium hypochlorite solution with mucous membrane, skin and eyes as well as the penetration into the periapical tissue.
The safety data sheet can be downloaded from www.legeartis.de or requested by E-Mail from sicherheitsdatenblaetter@legeartis.de.

Mode of application

After removing the closing cap attach or connect by turning a suitable disposable syringe onto the ESD-insert of the bottle neck and raise required volume. Subsequent to the positioning of a suitable cannula, rinse the root canal according to the individual rinse protocol.
For application of other rinsing methods, the solution can be taken via the opening for pouring out.
During the preparation of the root canal, the latter is rinsed with HISTOLITH NaOCl 1% by using suitable instruments (e.g. syringe with irrigation cannula, ultrasonic-activated or vibration-activated rinsing device) and with a qualified method (slow application without pressure, removing the rinsing solution by suction, protection of gingiva and oral mucosa by use of rubber dam) each time a new size of instrument is used, until treatment is complete.

Further notice

Sodium hypochlorite has a tissue-dissolving effect in particular on necrotic tissue. In case of an inappropriate application (e.g. pressing the solution through the apex) this may cause damages on vital tissue.

If the root canal is previously rinsed with EDTA-solution for widening the root canal and subsequently purified with HISTOLITH NaOCl 1%, this combination can be used to remove the smear layer during the root canal preparation.

The direct sequential use (without intermediate rinsing) of chlorhexidine and NaOCl solutions should be avoided.

In order to prevent pain and swelling through residues of sodium hypochlorite in the root canal the last rinsing can be done with a physiological saline solution.

Do not pour back residues of the solution into the bottle.

Keep out of reach of children.

Shelf life

Store in upright position, only in the original container.

HISTOLITH NaOCl 1% 200 ml has a shelf life of 2 years (unopened, storage at 2 – max. 25 °C).

HISTOLITH NaOCl 1% 500 ml has a shelf life of 3 years (unopened, storage at 2 – max. 25 °C).

The period of use after first opening is 6 months.

HISTOLITH NaOCl 1% shall not be used after the expiry date.

Storage instructions

The container has to be kept tightly closed, as otherwise a decrease of the NaOCl concentration has to be expected.

Administrative form and package sizes

200 ml solution Item number 0032342

500 ml solution Item number 0032343


Date of revision

2023-02

Symbol of “Expiry date”

Symbol of “Batch number”

Symbol of “Pay attention to the directions for use”


Symbol of “Manufacturer”

Symbol of “Catalogue number”
Shows Item number of the manufacturer

Symbol of “Temperature limitation”

Symbol of “Medical Device”

Symbol of “Unique device identification”


1. Handling of the ESD-syringe filling system with Luer or Luer Lock syringe


Remove the cap

Connect the syringe

Withdraw the desired volume

Remove the syringe

Close the cap

2. Pouring out the solution (without syringe) is also possible.

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