



Expert information

REMOT



Manufacturer

lege artis Pharma GmbH + Co KG
P. O. Box 60, D-72132 Dettenhausen
Breitwasenring 1, D-72135 Dettenhausen
Telephone +49 (0) 71 57 / 56 45 - 0
Fax +49 (0) 71 57 / 56 45 50
E-Mail: info@legeartis.de
www.legeartis.de

Name

REMOT prophylaxis paste
For removal of dental plaque, prophylaxis of caries and final polish.
For dental use.

Composition

Zirconium(IV)-oxide, silica, cocamidopropyl betaine, macrogols, sorbitol, glycerol, chlorhexidine digluconate and peppermint oil

Indications

REMOT removes plaque and therefore prevents caries and periodontal disease. Also suitable for final polish.

Notice

The well-balanced combination of abrasive materials effects a careful removal of the dental plaque without damaging the enamel.
REMOT has a pleasant refreshing taste, cleans teeth and makes them shiny.

Mode of application

REMOT will be applied with rubber cup and prophylaxis brushes.

Shelf life

REMOT shall not be used after the expiry date.

Administrative form and package size

75 g prophylaxis paste

Date of revision

January 2010

Symbol of "Expiry date"



Symbol of "Batch number"



0026317/Stückzahl0110/Bestelldatum