



# FOKALDRY



## The Turbo-DRYer

Solvent mixture for turbo cleaning  
in the dental practice & laboratory



## **FOKALDRY**

Solvent mixture for turbo cleaning, degreasing and drying cavities, root canals, tooth surfaces as well as crowns/bridges.

- ✓ Cleans, degreases and dries in one single step
- ✓ Dries quickly and residue-free
- ✓ Reduces germs
- ✓ Optimizes the adhesive bond with subsequent luting materials
- ✓ Contains acetone, isopropyl alcohol, petroleum ether, and diisopropyl ether

*Dosage form and package sizes: 500 ml solution (480 ml in storage bottle, 20 ml in further usage bottle plus 1 spout and 3 droppers) or 80 ml solution*

**lege artis Pharma  
GmbH + Co. KG**

Breitwasenring 1  
D-72135 Dettenhausen

+49 7157 5645-0  
+49 7157 5645-50  
info@legeartis.de  
www.legeartis.de

Proud sponsor of:



ARBEITSGEMEINSCHAFT  
DEUTSCHER  
KUNSTVEREINE

