



 SENSITIVE GINGIVA



20 ml / Medical Device (MD)

*UNCOMPLICATED
MOLDING*

ORBAT sensitive

The sensitive prep cleaner

25 % aluminium sulfate solution for gingival retraction, for simple and precise impression-taking.

- ✓ Enables a clean display of the preparation margin during impression-taking
- ✓ For sensitive gingiva
- ✓ For optimal tissue management during the cementation of direct and indirect restorations
- ✓ Compatible with impression materials
- ✓ For light bleeding
- ✓ No teeth discoloration: clear and colourless solution for high aesthetic demands

 lege artis

 lege artis



Get our
other
products
as well!



The successful prosthodontic classics of lege artis

- ✓ ORBAT sensitive ✓ ORBAT forte
- ✓ TECTOR ✓ Cp-CAP

lege artis Pharma
GmbH + Co. KG
Breitwasenring 1
D-72135 Dettenhausen

 +49 7157 5645-0
 +49 7157 5645-50
 info@legeartis.de
 www.legeartis.de

Proud sponsor of:



ARBEITSGEMEINSCHAFT
DEUTSCHER
KUNSTVEREINE



Version 02/2023

ACCORDING TO THE RULES OF ART

ACCORDING TO THE RULES OF ART

ORBAT forte
20 ml
Eisen(III)sulfat-
Flüssig-Gel

lege artis

20 ml / Medical Device (MD)

UNCOMPLICATED MOLDING

ORBAT forte

The forte prep cleaner

20 % aluminium sulfate solution for gingival retraction, for simple and precise impression-takings.

- ✓ Enables a clean display of the preparation margin during impression-taking
- ✓ For an optimal tissue management during the cementation of direct/indirect restorations
- ✓ Stops capillary bleeding
- ✓ Stops bleeding during pulpotomy procedures
- ✓ Must be rinsed off with water before impression-taking

TECTOR
10 ml
Suspension

lege artis

20 ml / Medical Device (MD)

PROTECTS THE DENTAL PULP

TECTOR

The pulp ProTECTOR

Calcium hydroxide suspension for protecting the dental pulp (open dentine tubules).

- ✓ Soothes in cases of grinding traumata
- ✓ Serves as a space-saving underfilling for shallow cavities and stumps before conventional cementation
- ✓ Neutralizes residual acids of luting material
- ✓ Ensures accurate fitting of restorations afterwards

Cp-CAP
15 g
Pulver

lege artis

20 ml / Medical Device (MD)

INDIVIDUAL MIXING

Cp-CAP

The CAP of the pulp

Zinc oxide eugenol cement with calcium hydroxide for maintaining the vitality of endangered dental pulp and for attaching long-term temporaries (2 components for individual mixing).

- ✓ For indirect pulp capping and for cementing temporary restorations
- ✓ Induces the formation of secondary and tertiary dentine
- ✓ Fast curing and x-ray visible/radiopaque